

## THERMOFLECT® HEAT REFLECTIVE TECHNOLOGY




 **thermoflect**® Heat Reflective Technology

Thermoflect® technology is the simple, safe and effective pre-warming protocol that maintains normothermia and prevents hypothermia.

It has been trusted by clinicians nationwide for over 20 years to protect at risk patients from the serious adverse events associated with unplanned hypothermia. Learn more at [www.Thermoflect.com](http://www.Thermoflect.com).



**NOT FOR  
USE IN MRI**

 **encompass**®

Better care starts  
with safety and comfort.

# THERMOFLECT® HEAT REFLECTIVE TECHNOLOGY



## Thermoflect® Heat Reflective Vests LATEX FREE

CAT NO	DESCRIPTION	SIZE	QTY/CASE
5130-200	Heat Reflective Vest	Small	25 EA
5130-300	Heat Reflective Vest	Medium	25 EA
5130-400	Heat Reflective Vest	Large	25 EA
5130-500	Heat Reflective Vest	X-Large	25 EA

## Thermoflect® Heat Reflective Staff Jackets LATEX FREE

CAT NO	DESCRIPTION	SIZE	QTY/CASE
5140-200	Silver Side Out	Small	25 EA
5140-250	Blue Side Out	Small	25 EA
5140-300	Silver Side Out	Medium	25 EA
5140-350	Blue Side Out	Medium	25 EA
5140-400	Silver Side Out	Large	25 EA
5140-450	Blue Side Out	Large	25 EA
5140-500	Silver Side Out	X-Large	25 EA
5140-550	Blue Side Out	X-Large	25 EA



## Thermoflect® Trauma Products LATEX FREE

CAT NO	DESCRIPTION	SIZE	QTY/CASE
5170-200	The Transport Cocoon™	Adult	10 EA
5170-150	The Transport Cocoon™	Pediatric	10 EA
5100-250	Papoose	5 to 10 lbs	50 EA

# THERMOFLECT® HEAT REFLECTIVE TECHNOLOGY



## Thermoflect® Adult Patient Apparel LATEX FREE

CAT NO	DESCRIPTION	SIZE	QTY/CASE
5110-100	Cap, Bouffant	Adult	100 EA
5190-600	Patient Gown, Adult	Small	25 EA
5145-500	Patient Gown, Adult	Universal	25 EA
5120-350	Leg Coverings	34" Long	50 PR
5115-200	Footies	Universal	25 PR

## Thermoflect® Pediatric Patient Apparel LATEX FREE

CAT NO	DESCRIPTION	SIZE	QTY/CASE
5110-250	Cap With Stockinette Band	Neonatal	100 EA
5110-200	Cap With Stockinette Band	Infant	100 EA
5110-150	Cap, Bouffant	Youth	50 EA
5180-200	Jacket (2-3)	X-Small	50 EA
5180-300	Jacket (4-5)	Small	50 EA
5180-400	Jacket (6-8)	Medium	50 EA
5180-500	Jacket (10-12)	Large	50 EA
5185-200	Pants (2-3)	X-Small	50 EA
5185-300	Pants (4-5)	Small	50 EA
5185-400	Pants (6-8)	Medium	50 EA
5185-500	Pants (10-12)	Large	50 EA
5190-200	Pediatric Gown (2-4)	Small	25 EA
5190-300	Pediatric Gown (6-8)	Medium	25 EA
5190-500	Pediatric Gown (10-12)	Large	25 EA



## Thermoflect® Patient Care LATEX FREE

CAT NO	DESCRIPTION	SIZE	QTY/CASE
5100-200	Blanket, Small	36" x 40"	60 EA
5100-400	Blanket, Medium	48" x 48"	50 EA
5100-700	Blanket, Large	48" x 84"	25 EA
5200-400	Procedure Pak 400 Medium Blanket, Leg Coverings, Adult Cap	Universal	20 EA
5200-700	Procedure Pak 700 Large Blanket, Adult Cap	Universal	24 EA
5105-100	Pulse Oximeter Mitt	Universal	50 EA
5300-100	TherMuff™ Bag, Sterile	8.5" x 9"	50 EA



**Thermoflect Heat Reflective Technology is the simple, safe and effective solution for preventing hypothermia and maintaining normothermia.**

#### Evidence-Based Practices for Thermoflect Heat Reflective Technology

*The Effect Of Preoperative Reflective Hats And Jackets, And Intraoperative Reflective Blankets On Perioperative Temperature, The Internet Journal of Anesthesiology. 2003. Volume 6, Yiwei Sheng, Frank Zavisca, Elizabeth Schonlau, Renee Desmarattes, Edwin Herron, Randall Cork, Department of Anesthesiology, LSUHSC Shreveport*

*Assessing the Use of a Reflective Blanket in Managing Post Operative Hypothermia, Rush University Medical Center and Rush University College of Nursing, Chicago IL, Kimberly Humbarger BSN RN, Isabel Pecson RN, Ruth Kleinpell PhD RN FAAN FCCM*

*The Usefulness of Preoperative Warming, B. van Beek, CRNA, Anaesthesia Department - Queen Beatrix Hospital, Winterswijk, The Netherlands, TSCI Abstract, October 2010*

*Thermal Efficiency of Pre-warmed Cotton, Reflective, and Forced Warm Air Inflatable Blankets in Trauma Patients, Emergency Nurses Association, 2000 Conference, Research Paper Presentation*

*Maintenance of Normothermia During Relatively Minor Orthopedic Surgical Procedures On the Extremities, Leo Leliveld, CRNA, Best Practices Lit. No. AC11187, January 2011*

*A Comparative Study of The Thermal Efficiency Of Thermo-Lite Clothing Versus Pre-Warmed Cotton Blankets In Preventing Intraoperative Hypothermia, By Joanne A. Kelley R.N., B.S.*



Better care starts with safety and comfort.

**800-284-4540 | [www.encompassgroup.net](http://www.encompassgroup.net)**

© 2011 Encompass Group, LLC All Rights Reserved Lit. No. AC10145 rev. 09/11